



**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for individuals
12 Years of Age and Older: Gray Cap and Label Border – Do NOT dilute**


Vial and carton labels with English-only labelling



VIAL (Gray Cap and Label Border)



Packaging no.	Pfizer-BioNTech COVID-19 Vaccine	NDC
LOT/MFG:	DO NOT DILUTE Age 12y & older Vial contains 6 doses of 0.3 mL For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. After first use store at 2 to 25°C (35 to 77°F) and discard after 6 hours. First use date and time:	

CARTON (10-Vial)



CARTON TEXT	
FRONT PANEL	
<p>Pfizer-BioNTech COVID-19 Vaccine <i>Suspension for Intramuscular Injection</i> DO NOT DILUTE 10 Multiple Dose Vials Each vial contains 6 doses of 0.3 mL For 12 years of age and older BIONTECH Pfizer</p> <p style="text-align: center;">TEAR HERE</p> 	
TOP PANEL	
<p>NDC Pfizer-BioNTech COVID-19 Vaccine <i>Suspension for Intramuscular Injection</i> DO NOT DILUTE 10 Multiple Dose Vials Each vial contains 6 doses of 0.3 mL For 12 years of age and older BIONTECH Pfizer Rx only</p>	
BACK PANEL	
<p>For storage and expiry information, see FDA-authorized Fact Sheet or scan QR code. Store in this carton to protect from light. After first use, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard 6 hours after first use. Contains no preservative. For use under Emergency Use Authorization</p> 	
BOTTOM PANEL	
<p>Manufactured for BioNTech Manufacturing GmbH An der Goldgrube 12 55131 Mainz, Germany Manufactured by Pfizer Inc. New York, NY 10017</p>	

SIDE PANEL

GTIN
LOT/MFG:

SIDE PANEL

Pfizer-BioNTech
COVID-19 Vaccine
Suspension for
Intramuscular Injection


DO NOT DILUTE

10 Multiple Dose Vials

Each vial contains

6 doses of 0.3 mL


For 12 years of age and older



**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for individuals
12 Years of Age and Older: Purple Cap and Label Border - DILUTE BEFORE USE**

Vial and carton labels with English-only labelling




VIAL (Purple Cap and Label Border)



Packaging no. Lot/Exp	Pfizer-BioNTech COVID-19 Vaccine After dilution, vial contains 6 doses of 0.3 mL For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. DILUTE BEFORE USE. Discard 6 hours after dilution when stored at 2 to 25°C (35 to 77°F). Dilution date and time:	NDC
--------------------------	---	-----

CARTON (195-Vial)

NO COPY

<p>NDC 59267-1000-2</p> <p>Pfizer-BioNTech COVID-19 Vaccine</p> <p><i>Suspension for Intramuscular Injection</i></p> <p>195 Multiple Dose Vials</p> <p>(after dilution each vial contains 6 doses of 0.3 mL)</p> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="text-align: center;">  </div> <div style="border: 1px solid red; padding: 2px; font-size: small;"> <p style="color: red;">FPO: QR Code</p> <p style="color: red;">https://www.pfi.sr/pfebntcovidvax</p> </div> <div style="text-align: center;">  </div> </div> <p style="font-size: x-small;"> STORAGE: Prior to dilution, store at -80°C to -60°C (-112°F to -76°F). Store in this carton to protect from light. DOSAGE AND ADMINISTRATION: After dilution, each vial contains 6 doses of 0.3 mL. See FDA-authorized Fact Sheet or scan QR code for information. MUST BE DILUTED BEFORE USE with sterile 0.9% Sodium Chloride Injection, USP (not supplied). After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard after 6 hours. Contains no preservative. For use under Emergency Use Authorization. Rx only </p>		<p style="font-size: x-small;"> Manufactured for BioNTech Manufacturing GmbH An der Goldgrube 12 55131 Mainz, Germany Manufactured by Pfizer Inc. New York, NY 10017 </p>	<p style="font-size: x-small; color: red;">359267100023</p>  <p style="font-size: x-small; color: red;">FPO UPC @ 80%</p>
NO COPY		<p style="font-size: large;">PAA173697</p>	
NO COPY		<p style="font-size: large;">PAA173697</p>	
<p>Pfizer-BioNTech COVID-19 Vaccine</p> <p><i>Suspension for Intramuscular Injection</i></p> <p>195 Multiple Dose Vials</p> <p>(after dilution each vial contains 6 doses of 0.3 mL)</p> <p style="font-size: x-small;"> STORAGE: Prior to dilution, store at -80°C to -60°C (-112°F to -76°F). Store in this carton to protect from light. DOSAGE AND ADMINISTRATION: After dilution, each vial contains 6 doses of 0.3 mL. See FDA-authorized Fact Sheet or scan QR code for information. </p>		<p>MUST BE DILUTED BEFORE USE with sterile 0.9% Sodium Chloride Injection, USP (not supplied). After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard after 6 hours. Contains no preservative.</p>	<p style="font-size: large;">PAA173697</p>
NO COPY		<p style="font-size: large;">PAA173697</p>	

CARTON TEXT

MAIN TOP PANEL

NDC 59267-1000-2

Pfizer-BioNTech COVID-19 Vaccine

Suspension for Intramuscular Injection

195 Multiple Dose Vials

(after dilution each vial contains 6 doses of 0.3 mL)

BIONTECH Pfizer



Manufactured for
BioNTech
Manufacturing GmbH
An der Goldgrube 12
55131 Mainz, Germany

Manufactured by
Pfizer Inc.
New York, NY 10017

STORAGE:

Prior to dilution, store at
-80°C to -60°C (-112°F to -76°F).
Store in this carton to protect from light.

Packaging no.

DOSAGE AND ADMINISTRATION:

After dilution, each vial contains 6 doses of 0.3 mL.
See FDA-authorized Fact Sheet or scan QR code for information.
MUST BE DILUTED BEFORE USE with sterile 0.9% Sodium
Chloride Injection, USP (not supplied).
After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).
Discard after 6 hours.
Contains no preservative.
For use under Emergency Use Authorization. Rx only

LOT:
EXP:

FRONT SIDE PANEL

Pfizer-BioNTech COVID-19 Vaccine

Suspension for Intramuscular Injection

195 Multiple Dose Vials

(after dilution each vial contains 6 doses of 0.3 mL)

STORAGE: Prior to dilution, store at
-80°C to -60°C (-112°F to -76°F).
Store in this carton to protect from light.

DOSAGE AND ADMINISTRATION:

After dilution, each vial contains 6 doses of 0.3 mL.
See FDA-authorized Fact Sheet or scan QR code for information.

MUST BE DILUTED BEFORE USE with sterile
0.9% Sodium Chloride Injection, USP (not supplied).
After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).
Discard after 6 hours.
Contains no preservative.


Packaging no.

LOT:
EXP:

**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for children
5 Years to <12 Years: **Orange Cap and Label Border** - **DILUTE PRIOR TO USE****

Vial and carton labels with English-only labelling

VIAL (Orange Cap and Label Border)



Packaging no.	Pfizer-BioNTech COVID-19 Vaccine	NDC:
LOT/MGF:	DILUTE PRIOR TO USE Age 5y to <12y After dilution – 10 doses of 0.2 mL For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. After dilution store at 2 to 25°C (35 to 77°F) and discard after 6 hours. Dilution date and time:	

CARTON (10-Vial)



CARTON TEXT

FRONT PANEL

Pfizer-BioNTech COVID-19 Vaccine

Suspension for Intramuscular Injection

DILUTE PRIOR TO USE


10 Multiple Dose Vials

After dilution each vial contains

10 doses of 0.2 mL

For age **5** years to **<12** years

BIONTECH Pfizer



TOP PANEL

NDC 59267-1055-4

**Pfizer-BioNTech
COVID-19 Vaccine**

Suspension for Intramuscular Injection

DILUTE PRIOR TO USE

10 Multiple Dose Vials

After dilution each vial contains 10 doses of 0.2 mL

For age **5** years to **<12** years

BIONTECH Pfizer Rx only

BACK PANEL

For storage and expiry information, see FDA-authorized Fact Sheet or scan QR code. Store in this carton to protect from light.

MUST BE DILUTED BEFORE USE with sterile

0.9% Sodium Chloride Injection, USP (not supplied).

After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).

Discard 6 hours after dilution.

Contains no preservative.

For use under Emergency Use Authorization



BOTTOM PANEL

Manufactured for
BioNTech
Manufacturing GmbH
An der Goldgrube 12
55131 Mainz, Germany
Manufactured by
Pfizer Inc.
New York, NY 10017

SIDE PANEL

GTIN
LOT/MFG:

SIDE PANEL

**Pfizer-BioNTech
COVID-19 Vaccine**

*Suspension for
Intramuscular Injection*

DILUTE PRIOR TO USE

10 Multiple Dose Vials

After dilution each

vial contains

10 doses of **0.2 mL**

For age **5** years to **<12** years

