



Vial and carton labels for COMIRNATY Original & Omicron BA.4/BA.5 (COVID-19 mRNA Vaccine, Bivalent [Original and Omicron BA.4/BA.5]) for individuals 12 years of age and older with English-only labelling

Vial Label



* Depending on product supply chain, the label component code number may be different.

Carton Label (10-vials)

Pfizer-BioNTech
COVID-19 Vaccine, Bivalent
 Original and Omicron BA.4/BA.5
 Suspension for Intramuscular Injection
DO NOT DILUTE
10 Multiple Dose Vials
 Each vial contains 6 doses of 0.3 mL
 For 12 years of age and older
 For use under Emergency Use Authorization.
 BIONTECH Pfizer

Pfizer-BioNTech
COVID-19 Vaccine, Bivalent
 Original and Omicron BA.4/BA.5
 Suspension for Intramuscular Injection
DO NOT DILUTE
10 Multiple Dose Vials
 Each vial contains 6 doses of 0.3 mL
 For 12 years of age and older

For storage and expiry information, see
 FDA-authorized Fact Sheet or scan QR code.
 Store in this carton to protect from light.
 After first use, store the vaccine at
 2°C to 25°C (35°F to 77°F).
 Discard 12 hours after first use.
 Contains no preservative.
 For use under Emergency Use Authorization.

Manufactured for
 BionTech
 Manufactured by
 An der Goldgrube 12
 55131 Mainz, Germany
 Pfizer Inc.
 New York, NY 10017
 MADE IN BELGIUM

359267 030429
 PAA194853

GTIN: 00989837030429
 LOT/MFG:

TEAR HERE

Rx only

NDC 59267-0304-2

FRONT PANEL:

Pfizer-BioNTech
COVID-19 Vaccine, Bivalent
Original and Omicron BA.4/BA.5
Suspension for Intramuscular Injection

DO NOT DILUTE

10 Multiple Dose Vials

Each vial contains **6** doses of **0.3 mL**

For **12 years of age and older**

BIONTECH Pfizer



TEAR HERE*



TOP PANEL:

Pfizer-BioNTech
COVID-19 Vaccine, Bivalent
Original and Omicron BA.4/BA.5
Suspension for Intramuscular Injection

NDC 59267-0304-2

DO NOT DILUTE

10 Multiple Dose Vials

Each vial contains **6** doses of **0.3 mL**

For **12 years of age and older**

For use under Emergency Use Authorization.

BIONTECH

Pfizer

Rx only

BACK PANEL:

For storage and expiry information, see
FDA-authorized Fact Sheet or scan QR code.
Store in this carton to protect from light.

After first use, store the vaccine at
2°C to 25°C (35°F to 77°F).

Discard 12 hours after first use.

Contains no preservative.

For use under Emergency Use Authorization.



BOTTOM PANEL:


Manufactured for
BioNTech
Manufacturing GmbH
An der Goldgrube 12
55131 Mainz, Germany

Manufactured by
Pfizer Inc.
New York, NY 10017


[Barcode]
359267030429*

MADE IN BELGIUM*

PAA194853*

SIDE PANEL
GTIN LOT/MFG:
SIDE PANEL
<p>Pfizer-BioNTech COVID-19 Vaccine, Bivalent Original and Omicron BA.4/BA.5 <i>Suspension for</i> <i>Intramuscular Injection</i> DO NOT DILUTE 10 Multiple Dose Vials Each vial contains  6 doses of 0.3 mL</p> <p>For 12 years of age and older</p>

* Depending on product supply chain and resulting differences in carton configuration, there may be the following slight differences in text and text location on the carton:

- The text "MADE IN BELGIUM" may not appear on all cartons.
- The label component code and bar code numbers may be different.
- The carton may not have a tearable flap; thus the text "↙ TEAR HERE ↘" may not be present.
- The upward arrows () designating the carton storage orientation may appear on a different panel.