


**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for Children Aged 6 Months to Less than 5 Years: MAROON Cap and Label Border - DILUTE PRIOR TO USE**


**Vial and carton labels with English-only labelling**

**VIAL (Maroon Cap and Label Border)**



Packaging no. LOT/MFG:	<p><b>Pfizer-BioNTech COVID-19 Vaccine</b>  <b>DILUTE PRIOR TO USE Age 6m to &lt;5y</b>                  After dilution – <b>10</b> doses of <b>0.2 mL</b>                  For intramuscular use. Contains no preservative.                  For use under Emergency Use Authorization.                  After dilution store at 2 to 25°C (35 to 77°F) and                  discard after 12 hours.                  Dilution date and time:</p>	NDC:
---------------------------	--	------

**CARTON (10-Vial)**



**CARTON TEXT**

**FRONT PANEL**

**Pfizer-BioNTech COVID-19 Vaccine**  
*Suspension for Intramuscular Injection*


**DILUTE PRIOR TO USE**  
**10 Multiple Dose Vials**

After dilution each vial contains  
**10 doses of 0.2 mL**  
For age **6 months to <5 years**

BIONTECH Pfizer



**TEAR HERE**



**TOP PANEL**

**Pfizer-BioNTech**  
**COVID-19 Vaccine**  
*Suspension for Intramuscular Injection*

NDC 59267-0078-4

**DILUTE PRIOR TO USE**  
**10 Multiple Dose Vials**

After dilution each vial contains  
**10 doses of 0.2 mL**  
For age **6 months to <5 years**

BIONTECH Pfizer

Rx only

**BACK PANEL**

For storage and expiry information, see FDA-authorized Fact Sheet or scan QR code. Store in this carton to protect from light.

**MUST BE DILUTED BEFORE USE** with sterile  
0.9% Sodium Chloride Injection,  
USP (not supplied).

After dilution, store the vaccine at  
2°C to 25°C (35°F to 77°F).

Discard 12 hours after dilution.

Contains no preservative.

For use under Emergency Use Authorization



**BOTTOM PANEL**

Manufactured for  
BioNTech  
Manufacturing GmbH  
An der Goldgrube 12  
55131 Mainz, Germany

Manufactured by  
Pfizer Inc.  
New York, NY 10017

Packaging no.

MADE IN GERMANY

**SIDE PANEL**

GTIN  
LOT/MFG:

**SIDE PANEL**


**Pfizer-BioNTech**  
**COVID-19 Vaccine**  
*Suspension for*  
*Intramuscular Injection*

**DILUTE PRIOR TO USE**  
**10 Multiple Dose Vials**

After dilution each  
vial contains

**10** doses of **0.2 mL**


For age **6** months to **<5** years



**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for individuals  
12 Years of Age and Older: Gray Cap and Label Border – Do NOT dilute**

**Vial and carton labels with English-only labelling**

**VIAL (Gray Cap and Label Border)**



Packaging no.	<b>Pfizer-BioNTech COVID-19 Vaccine</b>	NDC
LOT/MFG:	<b>DO NOT DILUTE      Age 12y &amp; older</b> Vial contains <b>6</b> doses of <b>0.3 mL</b> For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. After first use store at 2 to 25°C (35 to 77°F) and discard after 6 hours. First use date and time:	

**CARTON (10-Vial)**



<b>CARTON TEXT</b>	
<b>FRONT PANEL</b>	
<p><b>Pfizer-BioNTech COVID-19 Vaccine</b>  <i>Suspension for Intramuscular Injection</i>  <b>DO NOT DILUTE</b>  <b>10 Multiple Dose Vials</b>  Each vial contains <b>6</b> doses of <b>0.3 mL</b>  <b>For 12 years of age and older</b>  BIONTECH                  Pfizer</p> <p style="text-align: center;">TEAR HERE</p> 	
<b>TOP PANEL</b>	
<p>NDC  <b>Pfizer-BioNTech COVID-19 Vaccine</b>  <i>Suspension for Intramuscular Injection</i>  <b>DO NOT DILUTE</b>  <b>10 Multiple Dose Vials</b>  Each vial contains  6 doses of 0.3 mL  <b>For 12 years of age and older</b>  BIONTECH                  Pfizer                                  Rx only</p>	
<b>BACK PANEL</b>	
<p>For storage and expiry information, see  FDA-authorized Fact Sheet or scan QR code.  Store in this carton to protect from light.  After first use, store the vaccine at  2°C to 25°C (35°F to 77°F).  Discard 6 hours after first use.  Contains no preservative.  For use under Emergency Use Authorization</p> 	
<b>BOTTOM PANEL</b>	
<p>Manufactured for  BioNTech  Manufacturing GmbH  An der Goldgrube 12  55131 Mainz, Germany  Manufactured by  Pfizer Inc.  New York, NY 10017</p>	

**SIDE PANEL**

GTIN  
LOT/MFG:

**SIDE PANEL**

**Pfizer-BioNTech**  
**COVID-19 Vaccine**  
*Suspension for*  
*Intramuscular Injection*


**DO NOT DILUTE**

**10 Multiple Dose Vials**

Each vial contains

6 doses of 0.3 mL


**For 12 years of age and older**



**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for individuals  
12 Years of Age and Older: Purple Cap and Label Border - DILUTE BEFORE USE**

**Vial and carton labels with English-only labelling**

**VIAL (Purple Cap and Label Border)**



Packaging no. Lot/Exp	<b>Pfizer-BioNTech COVID-19 Vaccine</b> <b>After dilution, vial contains 6 doses of 0.3 mL</b> For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. DILUTE BEFORE USE. Discard 6 hours after dilution when stored at 2 to 25°C (35 to 77°F). Dilution date and time:	NDC
--------------------------	---	-----

**CARTON (195-Vial)**

NO COPY

<p>NDC 59267-1000-2</p> <p><b>Pfizer-BioNTech COVID-19 Vaccine</b></p> <p><i>Suspension for Intramuscular Injection</i></p> <p><b>195 Multiple Dose Vials</b></p> <p>(after dilution each vial contains 6 doses of 0.3 mL)</p> <p><b>BIONTECH</b> </p> <p><b>STORAGE:</b> Prior to dilution, store at -80°C to -60°C (-112°F to -76°F). Store in this carton to protect from light.</p> <p><b>DOSAGE AND ADMINISTRATION:</b> After dilution, each vial contains 6 doses of 0.3 mL. See FDA-authorized Fact Sheet or scan QR code for information.</p> <p><b>MUST BE DILUTED BEFORE USE</b> with sterile 0.9% Sodium Chloride Injection, USP (not supplied). After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard after 6 hours. Contains no preservative.</p> <p>For use under Emergency Use Authorization. <span style="float: right;"><b>Rx only</b></span></p>		<p>Manufactured for BioNTech Manufacturing GmbH An der Goldgrube 12 55131 Mainz, Germany Manufactured by Pfizer Inc. New York, NY 10017</p>	<p style="writing-mode: vertical-rl; transform: rotate(180deg);">359267100023</p> <p style="writing-mode: vertical-rl; transform: rotate(180deg);">FPO UPC @ 80%</p>
<p style="color: red;">FPO: QR Code</p> <p style="color: red;"><a href="https://www.pfi.sr/pfebntcovidvax">https://www.pfi.sr/pfebntcovidvax</a></p>		<p><b>PAA173697</b></p>	
<p style="color: red;">NO COPY</p>		<p style="color: red;">NO COPY</p>	
<p><b>Pfizer-BioNTech COVID-19 Vaccine</b></p> <p><i>Suspension for Intramuscular Injection</i></p> <p><b>195 Multiple Dose Vials</b></p> <p>(after dilution each vial contains 6 doses of 0.3 mL)</p> <p><b>STORAGE:</b> Prior to dilution, store at -80°C to -60°C (-112°F to -76°F). Store in this carton to protect from light.</p> <p><b>DOSAGE AND ADMINISTRATION:</b> After dilution, each vial contains 6 doses of 0.3 mL. See FDA-authorized Fact Sheet or scan QR code for information.</p>		<p><b>MUST BE DILUTED BEFORE USE</b> with sterile 0.9% Sodium Chloride Injection, USP (not supplied). After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard after 6 hours. Contains no preservative.</p>	<p><b>PAA173697</b></p>
<p><b>LOT:</b></p> <p><b>EXP:</b></p>		<p style="color: red;"><b>IMPRINT AREA</b> (reads this way) <b>NO VARNISH</b> <b>NO COPY, NO INK,</b></p>	
<p style="color: red;">NO COPY</p>		<p style="color: red;">NO COPY</p>	

**CARTON TEXT**

**MAIN TOP PANEL**

NDC 59267-1000-2

**Pfizer-BioNTech COVID-19 Vaccine**

*Suspension for Intramuscular Injection*

**195 Multiple Dose Vials**

(after dilution each vial contains 6 doses of 0.3 mL)

BIONTECH Pfizer



Manufactured for  
BioNTech  
Manufacturing GmbH  
An der Goldgrube 12  
55131 Mainz, Germany

Manufactured by  
Pfizer Inc.  
New York, NY 10017

**STORAGE:**

Prior to dilution, store at  
-80°C to -60°C (-112°F to -76°F).  
Store in this carton to protect from light.

Packaging no.

**DOSAGE AND ADMINISTRATION:**

After dilution, each vial contains 6 doses of 0.3 mL.  
See FDA-authorized Fact Sheet or scan QR code for information.  
**MUST BE DILUTED BEFORE USE** with sterile 0.9% Sodium  
Chloride Injection, USP (not supplied).  
After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).  
Discard after 6 hours.  
Contains no preservative.  
For use under Emergency Use Authorization. Rx only

LOT:  
EXP:

**FRONT SIDE PANEL**

**Pfizer-BioNTech COVID-19 Vaccine**

*Suspension for Intramuscular Injection*

**195 Multiple Dose Vials**

(after dilution each vial contains 6 doses of 0.3 mL)

**STORAGE:** Prior to dilution, store at

-80°C to -60°C (-112°F to -76°F).

Store in this carton to protect from light.

**DOSAGE AND ADMINISTRATION:**

After dilution, each vial contains 6 doses of 0.3 mL.  
See FDA-authorized Fact Sheet or scan QR code for information.

**MUST BE DILUTED BEFORE USE** with sterile  
0.9% Sodium Chloride Injection, USP (not supplied).  
After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).  
Discard after 6 hours.  
Contains no preservative.

Packaging no.

LOT:  
EXP:

**COMIRNATY (Pfizer-BioNTech COVID-19 Vaccine) for children  
5 Years to <12 Years: Orange Cap and Label Border - DILUTE PRIOR TO USE**

**Vial and carton labels with English-only labelling**

**VIAL (Orange Cap and Label Border)**



Packaging no.	<b>Pfizer-BioNTech COVID-19 Vaccine</b>	NDC:
LOT/MGF:	<b>DILUTE PRIOR TO USE Age 5y to &lt;12y</b> After dilution – 10 doses of 0.2 mL For intramuscular use. Contains no preservative. For use under Emergency Use Authorization. After dilution store at 2 to 25°C (35 to 77°F) and discard after 6 hours. Dilution date and time:	

**CARTON (10-Vial)**



**CARTON TEXT**

**FRONT PANEL**

**Pfizer-BioNTech COVID-19 Vaccine**

*Suspension for Intramuscular Injection*

**DILUTE PRIOR TO USE**

**10 Multiple Dose Vials**

After dilution each vial contains

**10** doses of **0.2 mL**

For age **5** years to **<12** years

BIONTECH Pfizer



**TOP PANEL**

NDC 59267-1055-4

**Pfizer-BioNTech  
COVID-19 Vaccine**

*Suspension for Intramuscular Injection*

**DILUTE PRIOR TO USE**

**10 Multiple Dose Vials**

After dilution each vial contains 10 doses of 0.2 mL

For age **5** years to **<12** years

BIONTECH Pfizer Rx only

**BACK PANEL**

For storage and expiry information, see FDA-authorized Fact Sheet or scan QR code. Store in this carton to protect from light.

**MUST BE DILUTED BEFORE USE** with sterile

0.9% Sodium Chloride Injection, USP (not supplied).

After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).

Discard 6 hours after dilution.

Contains no preservative.

For use under Emergency Use Authorization



**BOTTOM PANEL**

Manufactured for  
BioNTech  
Manufacturing GmbH  
An der Goldgrube 12  
55131 Mainz, Germany  
Manufactured by  
Pfizer Inc.  
New York, NY 10017

**SIDE PANEL**

GTIN  
LOT/MFG:

**SIDE PANEL**

**Pfizer-BioNTech  
COVID-19 Vaccine**

*Suspension for  
Intramuscular Injection*

**DILUTE PRIOR TO USE**

10 Multiple Dose Vials

After dilution each

vial contains

**10** doses of **0.2 mL**

For age **5** years to **<12** years

