

Outer Carton



Antigen Outer Label



Inner Label

2.5 mL DIN 02521326

FPO

LOT AAXXXXXX
EXP YYYY-MM-DD
412614

covifenz
COVID-19 Vaccine
(plant-based VLP,
recombinant, adjuvanted)

I.M. 3.75 mcg / dose

See leaflet
Voir le dépliant
COVIDVLP.com

Vial 1 of 2

Mix with adjuvant before use.
Before mixing, store at 2 to 8 °C.
After mixing, vial contains 10 doses
of 0.5 mL, must be used within
6 hours and stored at 20 to 30 °C.

**Mélanger avec l'adjuvant
avant l'utilisation.** Avant le
mélange, conserver entre 2 et
8 °C. Après le mélange, le flacon
contient 10 doses de 0,5 mL, doit
être utilisé dans les 6 heures et
entreposé entre 20 à 30 °C.

Mixing time / Heure de mélange :

medicago

Partie sans vernis

Dieline non-imprimable