



Updated vial and carton labels for Pfizer-BioNTech COVID-19 Vaccine with English-only labelling

Vial



Pfizer-BioNTech COVID-19 Vaccine
After dilution, vial contains 6 doses of 0.3 mL
For intramuscular use. Contains no preservative.
For use under Emergency Use Authorization.
DILUTE BEFORE USE. Discard 6 hours after
dilution when stored at 2 to 25°C (35 to 77°F).
Dilution date and time:

Carton



Main Panel

Pfizer-BioNTech COVID-19 Vaccine

Suspension for Intramuscular Injection

195 Multiple Dose Vials

(after dilution each vial contains 6 doses of 0.3 mL)

STORAGE:

Prior to dilution, store at -80°C to -60°C (-112°F to -76°F).
Store in this carton to protect from light.

DOSAGE AND ADMINISTRATION:

After dilution, each vial contains 6 doses of 0.3 mL
See FDA-authorized Fact Sheet or scan QR code for information.
MUST BE DILUTED BEFORE USE with sterile 0.9% Sodium Chloride Injection,
USP (not supplied).
After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F). Discard after 6 hours.

Contains no preservative.
For use under Emergency Use Authorization.



Manufactured by
Pfizer Inc
New York, NY 10017
Manufactured for
BioNTech
Manufacturing GmbH
An der Goldgrube 12
55131 Mainz, Germany

MADE IN GERMANY

Side Panel

Pfizer-BioNTech COVID-19 Vaccine

Suspension for Intramuscular Injection

195 Multiple Dose Vials

(after dilution each vial contains 6 doses of 0.3 mL)

STORAGE:

Prior to dilution, store at -80°C to -60°C (-112°F to -76°F)

Store in this carton to protect from light.

DOSAGE AND ADMINISTRATION:

After dilution, each vial contains 6 doses of 0.3 mL

See FDA-authorized Fact Sheet or scan QR code for information.

MUST BE DILUTED BEFORE USE with sterile

0.9% Sodium Chloride Injection, USP (not supplied).

After dilution, store the vaccine at 2°C to 25°C (35°F to 77°F).

Discard after 6 hours.

Contains no preservative.